

## Engle Serial Number Locations

### Dental Chair:

Use touch pad and lift chair to maximum height and recline the chair back fully back. Look under the chair seat on the left hand side. There will be a Block assembly midway back under the seat. You will see a sticker with the chair Serial number.



### Dental Delivery Systems: 360, 320, 310, Carts, Side & Wall Mount.

Serial number is located on the underside of the hand piece delivery system behind the tubing.



For further questions contact your local Engle Rep or Engle's Technical Support Dept.  
1-800-547-1906

[www.engledental.com](http://www.engledental.com)

## Assistant Instrumentation: 360 unit

Serial Number is located under the arm coming off the back of the chair.



## Dental Light:

Serial Number is located on the back of the light head

**Location as of November 2017**



**Serial Number Location prior to November 2017**

The Serial Number is located on the light transformer in the j-box at the front of the chair.



Halogen Light Transformer



LED Light Transformer

